


Product Datasheet

Vinculin antibody (orb106393)

Description	Rabbit polyclonal to Vinculin. Vinculin binding sites have been identified on other cytoskeletal pro...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...