


Product Datasheet

CDC25C antibody (orb106413)

Description	Rabbit polyclonal to CDC25C. a member of the MPI phosphatase family. Activates the cyclin dependent ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...