


Product Datasheet

SREBP-1 antibody (orb106425)

Description	Rabbit polyclonal to SREBP-1. This gene encodes a transcription factor that binds to the sterol regu...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of JK whole cell l...