


Product Datasheet

p21 Cip1 antibody (orb106430)

Description	Rabbit polyclonal to p21 Cip1. This gene encodes a potent cyclin-dependent kinase inhibitor. The enc...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...