


Product Datasheet

GluR1 antibody (orb106440)

Description	Rabbit polyclonal to GluR1. GluR1 an integral membrane protein belonging to the glutamate-gated ion ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...