


Product Datasheet

S6 Ribosomal Protein antibody (orb106461)

Description	Rabbit polyclonal to S6 Ribosomal Protein. Ribosomes, the organelles that catalyze protein synthesis...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...