


Product Datasheet

PkC theta antibody (orb106481)

Description	Rabbit polyclonal to PkC theta. Protein kinase C (PKC) is a family of serine- and threonine-specific...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...