

Product Datasheet

CPI17alpha antibody (orb106518)

Description	Rabbit polyclonal to CPI17alpha. PPP1R14A a phosphorylation-dependent inhibitory subunit of protein ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.

