


Product Datasheet

SirT1 antibody (orb106519)

Description	Rabbit polyclonal to SirT1. This gene encodes a member of the sirtuin family of proteins, homologs t...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...