


Product Datasheet

GAPDH antibody (orb106555)

Description	Rabbit polyclonal to GAPDH. Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa