


Product Datasheet

Human CCL19/MIP-3 Beta ELISA Kit (orb106604)

Description	Human CCL19/MIP-3 Beta ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G22-S98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10993 MW
Uniprot ID	Q6IBD6
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 40-2000 pg/ml.
Expiration Date	Please enquire.



Human
CCL19/MIP-
3 beta
PicoKine
ELISA Ki...