


Product Datasheet

Rat sL-Selectin ELISA Kit (orb106606)

Description	Rat sL-Selectin ELISA Kit
Reactivity	Rat
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W39-N332
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	42471 MW
Uniprot ID	Q63762
Sensitivity	300 pg/ml
Expiration Date	Please enquire.



Rat sL-
Selectin
PicoKine
ELISA Kit
stand...