


Product Datasheet

Mouse Angiopoietin-2 ELISA Kit (orb106607)

Description	Mouse Angiopoietin-2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: M1-F496
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	56576 MW
Uniprot ID	035608
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 10-20 ng/ml.
Expiration Date	Please enquire.



Mouse
Angiopoietin-
2 PicoKine
ELISA Kit ...