


Product Datasheet

Rat CCL4/MIP-1 Beta ELISA Kit (orb106608)

Description	Rat CCL4/MIP-1 Beta ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-N92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	10234 MW
Uniprot ID	P50230
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 30 pg/ml.
Expiration Date	Please enquire.



Rat
CCL4/MIP-1
beta
PicoKine
ELISA Kit s...