


Product Datasheet

Mouse CCL17/TARC ELISA Kit (orb106609)

Description	Mouse CCL17/TARC ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-P93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11539 MW
Uniprot ID	F6R5P4
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 16-147 pg/ml.
Expiration Date	Please enquire.



Mouse
CCL17/TARC
PicoKine
ELISA Kit
stan...