


Product Datasheet

Human CCL18/PARC ELISA Kit (orb106610)

Description	Human CCL18/PARC ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A21-A89
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	9849 MW
Uniprot ID	P55774
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum is around 40 ng/ml.
Expiration Date	Please enquire.



Human
CCL18/PARC
PicoKine
ELISA Kit
stan...