


Product Datasheet

Human B7-1/CD80 ELISA Kit (orb106611)

Description	Human B7-1/CD80 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V35-N242
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	33048 MW
Uniprot ID	P33681
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Human B7-1/CD80 PicoKine ELISA Kit stand...