


Product Datasheet

Rat CD86/B7-2 ELISA Kit (orb106612)

Description	Rat CD86/B7-2 ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V29-I250
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	17153 MW
Uniprot ID	Q62810
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 100 pg/ml.
Expiration Date	Please enquire.



Rat
CD86/B7-2
PicoKine
ELISA Kit
standar...