


Product Datasheet

Mouse CXCL4/PF4 ELISA Kit (orb106613)

Description	Mouse CXCL4/PF4 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V30-S105
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11243 MW
Uniprot ID	Q9Z126
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:200000, concentration in serum and plasma is 10-40 ug/ml.
Expiration Date	Please enquire.



Mouse
CXCL4/PF4
PicoKine
ELISA Kit
stand...