


Product Datasheet

Human CXCL13 / BLC / BCA1 ELISA Kit (orb106615)

Description	Human CXCL13 / BLC / BCA1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V23-R94
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum, plasma (heparin, EDTA) and saliva.
MW	12664 MW
Uniprot ID	043927
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 31-302 pg/ml.
Expiration Date	Please enquire.



Human
CXCL13/BLC
PicoKine
ELISA Kit
stan...