


Product Datasheet

PLK3 antibody (orb106866)

Description	Goat polyclonal antibody to PLK3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	polo-like kinase 3 (aa320-333)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell line HeLa (calculated MW of 71.6kDa according to NP_004064.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PDRLPISSCVTPD
Clonality	Polyclonal
MW	71.6
Entrez	1263
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using PLK3 antibody.



Immunohistochemical staining of Human Kidney tissue using PLK3 antibody.