


Product Datasheet

VCAM1 antibody (orb106867)

Description	Goat polyclonal antibody to VCAM1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	VCAM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 85kDa band observed in Human Uterus lysates (calculated MW of 81.3kDa according to NP_001069.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	TQIDSPLSGKVRSE
Clonality	Polyclonal
MW	81.3; 71.3; 74.3
Entrez	7412
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Uterus ly...