

Product Datasheet

SGCE antibody (orb106901)

Description	Rabbit polyclonal to SGCE. Component of the sarcoglycan complex, a subcomplex of the dystrophin-glyc...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SGCE
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49851
Uniprot ID	O43556
NCBI	NP_003910.1 , NP_001092870.1 , NP_001092871.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

