

Product Datasheet

PNPLA2 antibody (orb106908)

Description	Rabbit polyclonal to PNPLA2. Catalyzes the initial step in triglyceride hydrolysis in adipocyte and ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PNPLA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55316
Uniprot ID	Q96AD5
NCBI	NP_065109.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

