

Product Datasheet

SAPS2 antibody (orb106910)

Description	Rabbit polyclonal to SAPS2. Regulatory subunit of protein phosphatase 6 (PP6). May function as a sca...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPP6R2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	104942
Uniprot ID	O75170
NCBI	NP_001229829.1, NP_001229828.1, NP_055493.2, NP_001229827.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

