

Product Datasheet

OXGR1 antibody (orb106912)

Description	Rabbit polyclonal to OXGR1. Receptor for alpha-ketoglutarate. Seems to act exclusively through a G(q...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OXGR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38251
Uniprot ID	Q96P68
NCBI	NP_543008.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

