



Product Datasheet

ASB11 antibody (orb106916)

Description	Rabbit polyclonal to ASB11. May be a substrate-recognition component of a SCF-like ECS (Elongin-Cull...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ASB11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35367
Uniprot ID	Q8WXH4
NCBI	NP_001188512.1 , NP_001012428.1 , NP_543149.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis
of human
heart mouse
heart t...