

Product Datasheet

CTNNB1 antibody (orb106921)

Description Rabbit polyclonal to CTNNB1. Key downstream component of the canonical Wnt signaling pathway. In the...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target CTNNB1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 85497

Uniprot ID [P35222](#)

NCBI [NP_001895.1](#), [NP_001091679.1](#), [NP_001091680.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

