

Product Datasheet

ZMAT1 antibody (orb106922)

Description	Rabbit polyclonal to ZMAT1. Contains 3 matrin-type zinc fingers.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZMAT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74755
Uniprot ID	Q5H9K5
NCBI	NP_001269329.1 , NP_001269330.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

