

Product Datasheet

Mast1 antibody (orb106925)

Description	Rabbit polyclonal to Mast1. Appears to link the dystrophin/utrophin network with microtubule filamen...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Mast1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	170997
Uniprot ID	Q9R1L5
NCBI	NP_064329.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

