


Product Datasheet

GAPDH monoclonal antibody (orb10693)

Description	GAPDH monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Full length GAPDH protein purified from rabbit muscle
Target	GAPDH (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	38 kDa
Uniprot ID	P04406
Dilution Range	WB=1:1000-5000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Mouse brain lysate shows GAPDH expressio...