

Product Datasheet

TAS2R43 antibody (orb106931)

Description	Rabbit polyclonal to TAS2R43. Gustducin-coupled receptor implicated in the perception of bitter com...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TAS2R43
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35599
Uniprot ID	P59537
NCBI	NP_795365.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

