

Product Datasheet

HLA-DOA antibody (orb106934)

Description Rabbit polyclonal to HLA-DOA. Important modulator in the HLA class II restricted antigen presentatio...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target HLA-DOA

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 27599

Uniprot ID [P06340](#)

NCBI [NP_002110.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

