

Product Datasheet

p53 antibody (orb106938)

Description	Rabbit polyclonal antibody to p53
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TP53
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43653
Uniprot ID	P04637
NCBI	NP_001119586.1, NP_001119590.1, NP_001119585.1, NP_001263625.1, NP_001263626.1, NP_001119584.1, NP_001263624.1, NP_00126, NP_001119589.1, NP_001263627.1, NP_001119587.1, NP_000537.3, NP_001263628.1, NP_001119588.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

