

Product Datasheet

SGSM1 antibody (orb106945)

Description	Rabbit polyclonal to SGSM1. The function of this protein remains unknown.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SGSM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	129718
Uniprot ID	Q2NKQ1
NCBI	NP_001091967.1 , NP_597711.1 , NP_001091968.1 , NP_001035037.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

