

Product Datasheet

STK25 antibody (orb106946)

Description	Rabbit polyclonal to STK25. Oxidant stress-activated serine/threonine kinase that may play a role in...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	STK25
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48112
Uniprot ID	O00506
NCBI	NP_001269236.1 , NP_001258907.1 , NP_006365.2 , NP_001258906.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

