

Product Datasheet

PHKA2 antibody (orb106949)

Description	Rabbit polyclonal to PHKA2. Phosphorylase b kinase catalyzes the phosphorylation of serine in certai...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PHKA2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	138408
Uniprot ID	P46019
NCBI	NP_000283.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

