

Product Datasheet

RECQL antibody (orb106950)

Description	Rabbit polyclonal to RECQL. DNA helicase that may play a role in the repair of DNA that is damaged b...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RECQL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	73457
Uniprot ID	P46063
NCBI	NP_116559.1 , NP_002898.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

