

Product Datasheet

CPVL antibody (orb106952)

Description	Rabbit polyclonal to CPVL. May be involved in the digestion of phagocytosed particles in the lysosom...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CPVL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54164
Uniprot ID	Q9H3G5
NCBI	NP_061902.2 , NP_112601.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

