

Product Datasheet

GABRQ antibody (orb106974)

Description	Rabbit polyclonal to GABRQ. GABA, the major inhibitory neurotransmitter in the vertebrate brain, med...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GABRQ
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	71988
Uniprot ID	Q9UN88
NCBI	NP_061028.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

