

Product Datasheet

CWC27 antibody (orb106981)

Description	Rabbit polyclonal to CWC27. PPIases accelerate the folding of proteins (By similarity).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CWC27
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53847
Uniprot ID	Q6UX04
NCBI	NP_005860.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

