

Product Datasheet

ZNF12 antibody (orb106986)

Description	Rabbit polyclonal to ZNF12. Transcriptional repressor which suppresses activation protein 1 (AP-1)-a...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF12
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81202
Uniprot ID	P17014
NCBI	NP_008887.2 , NP_057349.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

