

Product Datasheet

CCDC155 antibody (orb106987)

Description	Rabbit polyclonal to CCDC155. The function of this protein remains unknown.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KASH5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62783
Uniprot ID	Q8N6LO
NCBI	NP_653289.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

