

Product Datasheet

TTC5 antibody (orb106989)

Description	Rabbit polyclonal to TTC5. Adapter protein involved in p53/TP53 response that acts by regulating and...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TTC5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48928
Uniprot ID	Q8N0Z6
NCBI	NP_612385.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

