

Product Datasheet

LRP2BP antibody (orb106991)

Description Rabbit polyclonal to LRP2BP. May act as an adapter that regulates LRP2 function.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target LRP2BP

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 39780

Uniprot ID [Q9P2M1](#)

NCBI [NP_060879.2](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

