

Product Datasheet

LRRC61 antibody (orb106992)

Description	Rabbit polyclonal to LRRC61. LRRC61 contains two Leucine-rich repeat; these are short sequence motif...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LRRC61
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28014
Uniprot ID	Q9BV99
NCBI	NP_001136400.1 , NP_076431.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

