

Product Datasheet

OXCT2 antibody (orb106994)

Description Rabbit polyclonal to OXCT2. Key enzyme for ketone body catabolism. Transfers the CoA moiety from suc...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target OXCT2

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 56140

Uniprot ID [Q9BYC2](#)

NCBI [NP_071403.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

