

Product Datasheet

SPAG6 antibody (orb106999)

Description	Rabbit polyclonal to SPAG6. Important for structural integrity of the central apparatus in the sperm...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SPAG6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55476
Uniprot ID	O75602
NCBI	NP_036575.1 , NP_001240784.1 , NP_758442.1 , NP_001240783.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

