

Product Datasheet

CLEC19A antibody (orb107001)

Description	Rabbit polyclonal to CLEC19A. The function of this protein remains unknown.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CLEC19A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15448
Uniprot ID	Q6UXS0
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

