

Product Datasheet

GINS3 antibody (orb107002)

Description	Rabbit polyclonal to GINS3. The GINS complex plays an essential role in the initiation of DNA replic...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GINS3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24535
Uniprot ID	Q9BRX5
NCBI	NP_001119602.1 , NP_073607.2 , NP_001119601.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

